




ExagonTM

Ceiling Pendant






The Project

The center of this project, is the need to place the patient and the room staff on the first place.

Exagon can accommodate medical supplies, monitors and ventilation machines. High load capacity ensures that vital equipment such as endoscopy and anesthesia machines are easily accessible. Customized workstations can be created easily and quickly, ensuring that everything is within comfortable reach.





Vertical

The guiding idea is the space saving of the vertical unit.



Modular

Exagon can accommodate medical supplies, monitors and ventilation machines. High load capacity ensures that vital equipment such as endoscopy and anesthesia machines are easily accessible.





• Lighting

The wall unit can be equipped with dimmable LEDs for ambient light, and with RGB units for better intraoperative contrast.



• Innovations

The absence of open surfaces and the integrated joints have allowed the filing of two international patents, testifying to the new line drawn by Exagon.



includes two new international patents



Exagon™
Ceiling Pendant

**“Exagon is what we were missing.
Finally more organization and efficiency in the operating room.”**



Materials

Its innovative project, using aerospace aluminum alloy, ensure high load capacity, up to 370 kgs.

Shelves



Standard shelves with drawer 50cmx50cm,
antibacterial paint,
led buttons for brakes,
double 25x10 inox steel eurorail



Medical glass shelves 50cmx50cm,
Antibacterial paint ,
Double 25x10 inox steel eurorail



Corian shelves 50cmx50cm,
Antibacterial paint ,
Double 25x10 inox steel eurorail

Data Sheet

Color: White or Customizable Color on Request

Paint: Antibacterial

Gas Outlets: All EU Standards Available

Arm Length **(X1)** : 800/ 900/ 1000/1200 mm

Maximum arm load **(X1)**: 371/ 320/ 279/ 216 Kg

Arm Length **(X2)** : 800/600/ 800 mm

Arm Length **(X1) + (X2)** : 1600/ 1600/ 2000 mm

Maximum arm load **(X1) + (X2)**: 135/ 135/ 83 Kg



INOX steel medical rails: Custom lenght on client need

Data sockets: HDMI, Rj45, Rj11, USB, Equipotential and varius data service available

Electrical sockets: all EU standards available

Lighting: LEDs for ambient light, RGB units for better intraopeartoire contrast.





Exagon™

Ceiling Pendant





Contacts

Mail to: info@dzmedicalesrl.com

Web: www.dzmedicalesrl.com


Microsite: www.exagon.dzmedicale.info



Exagon™

Ceiling Pendant

**“in moments of maximum concentration,
a space cleverly designed can save a life.”**



Exagon™
Ceiling Pendant



REV00 del
01/09/2025