


lineis™






lineis.

The bed head unit designed to bring a unique and innovative design to hospitals around the world.



Antibacterial Paint



lighting



Direct




Indirect



Night



Dimmer



LED



Electrical sockets,
all the standards available.

lineis, an innovative project, created by the union of only two high-strength aluminum profiles, an illuminated surface entirely integrated into the profile, which creates clean lines and harmonious shapes, allowing a more effective sanitation of the medical device.





Infrared remote control (not incl.)



The bed head unit can then be integrated with a 30x10 or 25x10 accessory bar, with electrical sockets on the lower part and switches on request. The design was also in agreement with the maintenance departments of the hospitals, to make it easy to inspect even after installation.



INOX steel medical rail:
30 x 10 // 25 x 10

data sheet

Weight: 6 Kg per meter

Color: White, Black + customizable colors on request

Paint: Antibacterial

Gas Outlets: Variables for type and number

Section Measurements: Length based on needs, or on customer request

INOX steel medical rail: 30 x 10 // 25 x 10

Closures: lateral without screws

GAS socket: DIN and special

Data sockets: HDMI, Rj45, Rj11, USB, Equipotential and various data service available

Electrical sockets: all the standards available

Lighting: Full-length LED (patient)
Double full-length LED light (ambient + reading)

Remote control: infrared (not included but can be ordered separately)

lineis



is pure design



REV00 del
01/09/2025

for more information visit the product microsite:

www.lineis.dzmedicale.info

